

Product Datasheet

Integrin alpha 1/CD49a Antibody (CL7207) - BSA Free NBP2-76478

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-76478

Updated 9/22/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-76478



NBP2-76478

Integrin alpha 1/CD49a Antibody (CL7207) - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL7207
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS, pH 7.2, containing 40% glycerol

Product Description	
Description	Novus Biologicals Mouse Integrin alpha 1/CD49a Antibody (CL7207) - BSA Free (NBP2-76478) is a monoclonal antibody validated for use in IHC and WB. Anti-Integrin alpha 1/CD49a Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3672
Gene Symbol	ITGA1
Species	Human
Reactivity Notes	Mouse (89%).
Immunogen	This antibody was developed using a synthetic peptide derived from P56199, with the exact immunogen sequence remaining proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Block/Neutralize, Knockdown Validated
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500, Knockdown Validated, Block/Neutralize Reported in scientific literature (PMID:31750298).
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Immunohistochemistry-Paraffin: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478] - Analysis in human prostate and skeletal muscle tissues. Corresponding ITGA1 RNA-seq data are presented for the same tissues.




Western Blot: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478]
 - Analysis in U-138MG cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-ITGA1 antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH.




Immunohistochemistry-Paraffin: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478] - Staining of human placenta shows strong membranous positivity in blood vessels and loose connective tissue of chorionic villi.




Immunohistochemistry-Paraffin: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478] - Staining of human prostate shows strong membranous positivity in smooth muscle cells.



Immunohistochemistry-Paraffin: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478] - Staining of human skeletal muscle shows no positivity in striated muscle fibers as expected.



Immunohistochemistry-Paraffin: Integrin alpha 1/CD49a Antibody (CL7207) [NBP2-76478] - Staining of human stomach shows moderate membranous positivity in glandular cells.



Publications

Wu D, Witt RL, Harrington DA, Farach-Carson MC Dynamic Assembly of Human Salivary Stem/Progenitor Microstructures Requires Coordinated alpha1 beta1 Integrin-Mediated Motility Front Cell Dev Biol 2019-10-16 [PMID: 31750298] (B/N, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-76478

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-76478

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

