

Product Datasheet

Polyethylene Glycol Antibody (RM105) NBP2-77399

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-77399

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-77399



NBP2-77399**Polyethylene Glycol Antibody (RM105)**

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	RM105
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	50% Glycerol/PBS, 1% BSA


Product Description	
Description	Novus Biologicals Rabbit Polyethylene Glycol Antibody (RM105) (NBP2-77399) is a recombinant monoclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Species	All Species
Specificity/Sensitivity	This antibody reacts to the methoxy group of Polyethylene glycol (PEG). It is specific for methoxypoly(ethylene glycol).
Immunogen	KLH conjugated Polyethylene Glycol (PEG) with terminal methoxy group

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry
Recommended Dilutions	Western Blot 0.05 ug/ml - 1 ug/ml, ELISA 0.01 ug/ml-0.3 ug/ml, Immunohistochemistry 0.5ug/ml - 2ug/ml, Immunohistochemistry-Paraffin




Images


Western Blot: Polyethylene Glycol Antibody (RM105) [NBP2-77399] - BSA and PEGylated BSA (mPEG 5 kDa) using 0.1mg/mL of antiPEG clone RM105.



Immunohistochemistry: Polyethylene Glycol Antibody (RM105) [NBP2-77399] - Mouse liver and kidney using 0.5mg/mL of anti-PEG Clone RM105. Mouse had been injected with PEG-BSA or BSA for 3 hours before sampling.



ELISA: Polyethylene Glycol Antibody (RM105) - Azide and BSA Free [NBP2-77399] - ELISA of three different PEGylated BSAs using Clone RM105, followed by an AP conjugated goat anti-rabbit IgG.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-77399

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77399

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

