

Product Datasheet

MafA Antibody (BLR067G) - BSA Free NBP2-78707

Unit Size: 100 uL

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-78707

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-78707



NBP2-78707


MafA Antibody (BLR067G) - BSA Free

Product Information	
Unit Size	100 uL
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 2 - 8 C / 1 year from date of receipt
Clonality	Monoclonal
Clone	BLR067G
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A or G purified from cell culture supernatant
Buffer	Borate Buffered Saline (BBS) pH 8.2.
Product Description	
Description	Novus Biologicals Rabbit MafA Antibody (BLR067G) - BSA Free (NBP2-78707) is a recombinant monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	389692
Gene Symbol	MAFA
Species	Human, Mouse
Immunogen	Immunogen was a peptide representing a region between residue 125 and 175 of human MAF bZIP transcription factor A using the numbering given in entry NP_963883.2 (Gene ID 389692).
Notes	This Product sold under License from Bethyl Laboratories, Inc.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100 - 1:500, Immunohistochemistry-Paraffin
Application Notes	Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. IHC and ICC/IF: Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections. Zinc-fixative (JB fix) is recommended to enhance staining.




Images


Western Blot: MafA Antibody (BLR067G) [NBP2-78707] - Detection of human and mouse MafA by WB of HEK293T transfected with myc tagged human or mouse MafA, MafB or c-Maf using this antibody, rabbit anti-MafB recombinant monoclonal [BLR046F], and rabbit anti-c-Maf recombinant monoclonal [BLR045F]. Secondary: HRP-conjugated goat anti-rabbit IgG (NB100-2332). Lower panel: Rabbit anti-Myc.




Immunocytochemistry/Immunofluorescence: MafA Antibody (BLR067G) [NBP2-78707] - Rabbit anti-MafA recombinant monoclonal antibody [BLR067G]. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.




Immunohistochemistry-Paraffin: MafA Antibody (BLR067G) [NBP2-78707] - Rabbit anti-MafA recombinant monoclonal antibody [BLR067G]. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.



Western Blot: MafA Antibody (BLR067G) [NBP2-78707] - Detection of mouse MafA by WB of Beta-TC-6, AlphaTC1 Clone 9, and NIH 3T3 lysate using this antibody. Secondary: HRP-conjugated goat anti-rabbit IgG (NB100-2332).



Immunohistochemistry-Paraffin: MafA Antibody (BLR067G) [NBP2-78707] - Detection of human MafA (red) in FFPE pancreatic carcinoma by IHC-IF using this antibody. Secondary: HRP-conjugated goat anti-rabbit IgG (NB100-2332). Substrate: Opal™. Counterstain: DAPI (blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-78707

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-78707

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

