

# Product Datasheet

## GAD2/GAD65 Antibody (GAD2/1960) NBP2-79804-100ug

Unit Size: 100 ug

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-79804](http://www.novusbio.com/NBP2-79804)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-79804](http://www.novusbio.com/reviews/destination/NBP2-79804)



**NBP2-79804-100ug**

GAD2/GAD65 Antibody (GAD2/1960)

| Product Information     |                          |
|-------------------------|--------------------------|
| Unit Size               | 100 ug                   |
| Concentration           | 0.2 mg/ml                |
| Storage                 | Store at 4C.             |
| Clonality               | Monoclonal               |
| Clone                   | GAD2/1960                |
| Preservative            | 0.05% Sodium Azide       |
| Isotype                 | IgG2b Kappa              |
| Purity                  | Protein A or G purified  |
| Buffer                  | 10 mM PBS with 0.05% BSA |
| Target Molecular Weight | 65 kDa                   |


| Product Description |  |
|---------------------|--|
| Description         | 200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP2-79938)<br><br>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C. |
| Host                | Mouse  |
| Gene ID             | 2572   |
| Gene Symbol         | GAD2   |
| Species             | Human  |
| Marker              | GABAergic Neuronal Marker  |
| Immunogen           | Recombinant human GAD2/GAD65 protein fragment (around aa 6-99) (exact sequence is proprietary) (Uniprot: Q05329)   |

| Product Application Details |   |
|-----------------------------|---|
| Applications                | Immunohistochemistry, Immunohistochemistry-Paraffin, Protein Array  |
| Recommended Dilutions       | Immunohistochemistry 0.1-0.2 ug/ml, Immunohistochemistry-Paraffin 0.1-0.2 ug/ml, Protein Array  |
| Application Notes           | Immunohistochemistry (Formalin-fixed): 0.1-0.2ug/ml for 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined. |




## Images

Immunohistochemistry-Paraffin: GAD2/GAD65 Antibody (GAD2/1960) [NBP2-79804] - Formalin-fixed, paraffin-embedded human Pancreas stained with GAD2/GAD65 Antibody (GAD2/1960).



Protein Array: GAD2/GAD65 Antibody (GAD2/1960) [NBP2-79804] - Analysis of Protein Array containing more than 19,000 full-length human proteins using GAD2/GAD65 Antibody (GAD2/1960) Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (Monoclonal Antibody) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProt(TM) array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt(TM) are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a Monoclonal Antibody to its intended target. A Monoclonal Antibody is considered to specific to its intended target, if the Monoclonal Antibody has an S-score of at least 2.5





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-79804-100ug**

---

|                    |   |
|--------------------|---|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-43317-0.5mg   | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)    |
| H00002572-P02-10ug | Recombinant Human GAD2/GAD65 GST (N-Term) Protein       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-79804](http://www.novusbio.com/reviews/submit/NBP2-79804)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

