

Product Datasheet

Human Connexin 31/GJB3 ELISA Kit (Colorimetric) NBP2-79813

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-79813

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-79813




NBP2-79813**Human Connexin 31/GJB3 ELISA Kit (Colorimetric)**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	2707
Gene Symbol	GJB3
Species	Human
Specificity/Sensitivity	This kit recognizes Human CX31 in samples. No Significant cross-reactivity or interference between Human CX31 and analogues was observed.
Kit Components	Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Manual
Standard Curve Range	0.31-20 ng/mL
Sensitivity	0.19 ng/mL
Inter Assay	CV% < 4.78%
Intra Assay	CV% < 5.38%
Assay Type	Sandwich-ELISA
Spiking Recovery	86-108%
Suitable Sample Type	Serum, plasma, and other biological fluids
Sample Volume	100 ul
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human Connexin 31/GJB3 ELISA Kit (Colorimetric) [NBP2-79813] - Samples were spiked with high concentrations of Human Connexin 31/GJB3 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human Connexin 31/GJB3 ELISA Kit (Colorimetric) [NBP2-79813] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-79813

H00002707-P01-10ug AF748	Recombinant Human Connexin 31/GJB3 GST (N-Term) Protein E-Cadherin Antibody [Unconjugated]
H00002707-M01	Connexin 31/GJB3 Antibody (3B4-1B3)
NBP1-85237	SLC26A4 Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-79813

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

