

Product Datasheet

PBEF/Visfatin/NAMPT Antibody (OMNI379) - Azide and BSA Free NBP2-80036-100ug

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80036

Updated 9/11/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80036




NBP2-80036-100ug

PBEF/Visfatin/NAMPT Antibody (OMNI379) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OMNI379
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein G purified
Buffer	0.2 um-filtered solution in PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	10135
Gene Symbol	NAMPT
Species	Human, Mouse, Rat
Immunogen	Recombinant human PBEF/Visfatin/NAMPT (visfatin/PBEF).
Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, CyTOF-ready
Recommended Dilutions	Western Blot 1:2000 - 1:5000, Flow Cytometry 1:1000, ELISA 1:2000 - 1:10000, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/Immunofluorescence 1:500 - 1:1000, Immunoprecipitation 1:500 - 1:1000, Immunohistochemistry-Paraffin 1:100 - 1:500, CyTOF-ready

Images


Western Blot: PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036]
 - Western blot analysis of 3T3-L1 cell lysate using anti-Nampt (Visfatin/PBEF) mAb (OMNI379) at 1:2000 dilution. 1. IBMX. + DEX. - treated 3T3-L1 cell lysate. 2. IBMX. + DEX. + Insulin - treated 3T3-L1 cell lysate.




Immunocytochemistry/Immunofluorescence: PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036] - Immunofluorescent confocal imaging of HEK293 cells using anti-Nampt (Visfatin-PBEF), mAb (OMNI379). Nampt (Visfatin/PBEF) (green staining) is localized in the cytoplasm and nucleus (400x magnification). Nuclei are stained with Hoechst (blue staining).




Immunohistochemistry-Paraffin: PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036] - Immunohistochemical staining of Nampt using anti-Nampt (OMNI379) in human kidney tissue (1:500 dilution). Immunoperoxidase staining of formalin-fixed and paraffin-embedded (FFPE) human kidney showing cytoplasmic staining (200x, brown color).




Immunohistochemistry-Paraffin: PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036] - Immunohistochemical staining of Nampt using anti-Nampt (Visfatin/PBEF), mAb (OMNI379) in formalin-fixed and paraffin-embedded (FFPE) human colon tissue (15 ug/mL).



PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036] - Western blot analysis using anti-Nampt (Visfatin-PBEF), mAb (OMNI379) at 1:2000 dilution. 1. Human Nampt (Visfatin/PBEF) (His-tagged). 2. Human Nampt (Visfatin/PBEF) (FLAG(R)-tagged). 3. Mouse Nampt (Visfatin/PBEF) (His-tagged). 4. Rat Nampt (Visfatin/PBEF) (His-tagged). 5. Human Nampt (Visfatin/PBEF) (His-tagged). 6. LPS-treated human peripheral blood leukocyte lysate.



Immunoprecipitation: PBEF/Visfatin/NAMPT Antibody (OMNI379) [NBP2-80036] - Immunoprecipitation of recombinant human Nampt (Visfatin/PBEF) proteins using anti-Nampt (Visfatin-PBEF), mAb (OMNI379). Recombinant human Nampt (Visfatin/PBEF) proteins at different concentrations were precipitated. The precipitated proteins were separated by SDS-PAGE, electroblotted and visualized by western blot with anti-human Nampt (Visfatin/PBEF) pAb.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80036-100ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBC1-22597	Recombinant Mouse PBEF/Visfatin/NAMPT His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80036

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

