

Product Datasheet

Human NK3R/TACR3/Neurokinin B Receptor ELISA Kit (Colorimetric) NBP2-80294

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80294

Updated 9/23/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80294



NBP2-80294**Human NK3R/TACR3/Neurokinin B Receptor ELISA Kit (Colorimetric)**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.

Product Description	
Gene ID	6870
Gene Symbol	TACR3
Species	Human
Specificity/Sensitivity	This kit recognizes Human TACR3 in samples. No significant cross-reactivity or interference between Human TACR3 and analogues was observed.
Kit Components	Micro ELISA Plate(Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100X), Concentrated HRP Conjugate (100X), Reference Standard and Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25X), Substrate Reagent, Stop Solution, Plate Sealer, Product Description, Certificate of Analysis
Standard Curve Range	0.156-10 ng/mL
Sensitivity	0.094 ng/mL
Inter Assay	CV% < 4.46%
Intra Assay	CV% < 5.44%
Assay Type	Sandwich-ELISA
Spiking Recovery	86-104%
Suitable Sample Type	Serum, plasma and other biological fluids
Sample Volume	100 ul

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human NK3R/TACR3/Neurokinin B Receptor ELISA Kit (Colorimetric) [NBP2-80294] - Standard Reference Curve





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80294

NB300-102PEP	NK3R/TACR3/Neurokinin B Receptor Antibody Blocking Peptide
247-ILB-005/CF	IL-15 [Unconjugated]
NB300-102	NK3R/TACR3/Neurokinin B Receptor Antibody
MAB139-100	NKG2D/CD314 Antibody (149810) [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80294

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

