

# Product Datasheet

## SLC30A5 Antibody NBP2-80523

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-80523](http://www.novusbio.com/NBP2-80523)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-80523](http://www.novusbio.com/reviews/destination/NBP2-80523)



**NBP2-80523****SLC30A5 Antibody**

| <b>Product Information</b>         |                                                                                                                                                                                      |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                               |
| <b>Concentration</b>               | This product is unpurified. The exact concentration of antibody is not quantifiable.                                                                                                 |
| <b>Storage</b>                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                               |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                           |
| <b>Preservative</b>                | No Preservative                                                                                                                                                                      |
| <b>Reconstitution Instructions</b> | Reconstitute in 0.1 ml of sterile water. Centrifuge to remove any insoluble material. Glycerol may be added (1:1) for additional stability.                                          |
| <b>Purity</b>                      | Unpurified                                                                                                                                                                           |
| <b>Buffer</b>                      | Lyophilized from whole serum                                                                                                                                                         |
| <b>Product Description</b>         |                                                                                                                                                                                      |
| <b>Description</b>                 | Novus Biologicals Rabbit SLC30A5 Antibody (NBP2-80523) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>                        | Rabbit                                                                                                                                                                               |
| <b>Gene ID</b>                     | 64924                                                                                                                                                                                |
| <b>Gene Symbol</b>                 | SLC30A5                                                                                                                                                                              |
| <b>Species</b>                     | Human, Mouse, Rat                                                                                                                                                                    |
| <b>Immunogen</b>                   | A synthetic peptide from internal cytoplasmic region of human SLC30A5 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in rat and mouse.        |
| <b>Product Application Details</b> |                                                                                                                                                                                      |
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                                                                                                    |
| <b>Recommended Dilutions</b>       | Western Blot 1:2000, Immunohistochemistry 1:2000, Immunohistochemistry-Paraffin 1:2000                                                                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-80523**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| H00064924-P01-10ug | Recombinant Human SLC30A5 GST (N-Term) Protein      |
| M6000B-1           | IL-6 [HRP]                                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-80523](http://www.novusbio.com/reviews/submit/NBP2-80523)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

