

Product Datasheet

alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free NBP2-80569

Unit Size: 0.25 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80569

Updated 9/11/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80569



NBP2-80569

alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free


Product Information	
Unit Size	0.25 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	YOL1/34
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	50 kDa

Product Description	
Host	Rat
Gene ID	7846
Gene Symbol	TUBA1A
Species	Human, Mouse, Rat, Drosophila, Fungi, Yeast
Reactivity Notes	Cross reacts with Arabidopsis, Saccharomyces, Platyzoa, Ashbya, and Naegleria. Based on sequence similarity, is expected to react with: Birds and Mammals.
Immunogen	This alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free was developed against yeast tubulin.


Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Radioimmunoassay
Recommended Dilutions	Western Blot 1:100-1:2000, Flow Cytometry, ELISA 10 ug/ml, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Radioimmunoassay 1:100-1:2000, Flow (Intracellular)

Images


Western Blot: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - Total protein from human HeLa and A431, mouse 3T3 and rat PC12 cell lines was separated on a 12% gel by SDS-PAGE, transferred to PVDF membrane and blocked in 5% non-fat milk in TBST. The membrane was probed with 2.0 ug/ml anti-alpha Tubulin in block buffer and detected with an anti-mouse HRP secondary antibody using West Pico PLUS chemiluminescence detection reagent. Alpha tubulin molecular weight: 50 kDa. Image from the standard format of this antibody.




Immunocytochemistry/Immunofluorescence: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - PC12 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X PBS + 0.05% Triton-X100. The cells were incubated with anti-Tubulin Antibody [YOL1/34] at 5 ug/ml overnight at 4C and detected with an anti-rat Dylight 488 (Green) at a 1:500 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective. Image from the standard format of this antibody.




Flow Cytometry: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - An intracellular stain was performed on HeLa cells with alpha Tubulin Antibody [YOL1/34] NB100-1639AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 2.5 ug/mL for 30 minutes at room temperature. Both antibodies were directly conjugated to Alexa Fluor 647. Image using the Alexa Fluor 647 format of this antibody.




Immunocytochemistry/Immunofluorescence: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - HeLa cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X PBS + 0.05% Triton-X100. The cells were incubated with anti-Tubulin Antibody [YOL1/34] at 5 ug/ml overnight at 4C and detected with an anti-rat Dylight 488 (Green) at a 1:500 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective. Image from the standard format of this antibody.



Immunocytochemistry/Immunofluorescence: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - Representative images from Cdc14 release assay. Cells are assayed for spindle morphology & Cdc14 release phenotype by indirect immunofluorescence with anti-tubulin or anti-HA, which recognizes the HA-tagged version of Cdc14 in all strains. Cells are visualized by nuclear staining with DAPI & by differential interference contrast (DIC) optics for comparison. No release refers to dense nucleolar staining of Cdc14-HA, full release refers to cells with uniform Cdc14-HA localization throughout the nucleus, while partial release refers to cells with both nuclear & residual nucleolar staining. Image collected & cropped by CiteAb from the following publication (<https://dx.plos.org/10.1371/journal.pgen.1000588>), licensed under a CC-BY license. Image from the standard format of this antibody.



Immunocytochemistry/Immunofluorescence: alpha Tubulin Antibody (YOL1/34) - Azide and BSA Free [NBP2-80569] - NIH3T3 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X PBS + 0.05% Triton-X100. The cells were incubated with anti-Tubulin Antibody [YOL1/34] at 5 ug/ml overnight at 4C and detected with an anti-rat Dylight 488 (Green) at a 1:500 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80569

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NBP1-75398	Goat anti-Rat IgG (H+L) Secondary Antibody (Pre-adsorbed)
NBP2-21947-0.1mg	Rat IgG2a Isotype Control (2A3)
H00007846-Q01-10ug	Recombinant Human alpha Tubulin GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80569

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

