

Product Datasheet

EBI3 Antibody (10J763) - Azide and BSA Free NBP2-80705

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80705

Updated 8/21/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80705



NBP2-80705


EBI3 Antibody (10J763) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	10J763
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein G purified
Buffer	PBS
Product Description	
Host	Rat
Gene ID	10148
Gene Symbol	EBI3
Species	Mouse
Immunogen	Full-length recombinant mouse EBI3 protein
Product Application Details	
Applications	Western Blot, Flow Cytometry, CyTOF-ready
Recommended Dilutions	Western Blot 0.1-3 ug/ml, Flow Cytometry 1-3ug/10 ⁶ cells, CyTOF-ready




Images


Western Blot: EBI3 Antibody (10J763) - Azide and BSA Free [NBP2-80705] - Western blot testing of mouse EBI3 monoclonal antibody at A) 0.1 ug/ml on recombinant protein and B) 3 ug/ml on recombinant mouse IL-35 protein Goat anti-rat Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test. Image from the standard format of this antibody.



Flow Cytometry: EBI3 Antibody (10J763) - Azide and BSA Free [NBP2-80705] - Figure A: Intracellular analysis of mouse IL-35 in transfected and untransfected HEK 293 cells (Brefeldin A treated, 5 hours) using mouse EBI3 antibody at 3 ug/10⁶ cells. Figure B: Intracellular analysis of mouse EBI3 in stimulated and unstimulated RAW cells (LPS, 50 ng/ml, overnight) using mouse EBI3 antibody at 1.5 ug/10⁶ cells. Image from the standard format of this antibody.



Western Blot: EBI3 Antibody (10J763) - Azide and BSA Free [NBP2-80705] - Western blot analysis of mouse thymus tissue using EBI3 antibody at a concentration of 1 ug/ml. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80705

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NBP1-75398	Goat anti-Rat IgG (H+L) Secondary Antibody (Pre-adsorbed)
NBP1-43321-0.5mg	Rat IgG2a Kappa Light Chain Isotype Control (R2a)
NBP2-76293-100ug	Recombinant Human EBI3 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80705

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

