

Product Datasheet

PINK1 Antibody (8E10.1D6) - Azide and BSA Free NBP2-80910

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80910

Updated 8/21/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80910



NBP2-80910

PINK1 Antibody (8E10.1D6) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8E10.1D6
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	62.7 kDa


Product Description	
Host	Mouse
Gene ID	65018
Gene Symbol	PINK1
Species	Human, Mouse, Rat
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 29486776).
Immunogen	PINK1 antibody was developed using a synthetic peptide made to the human PINK1 protein sequence (between residues 100-250). [Swiss-Prot: Q9BXM7]

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, SDS-Page, Knockdown Validated
Recommended Dilutions	Western Blot 2-4 ug/ml, Immunohistochemistry 5 ug/ml, Immunocytochemistry/ Immunofluorescence 20-50 ug/ml, SDS-Page reported in scientific literature (PMID 27553674), Knockdown Validated
Application Notes	Unprocessed PINK1 is 63 kDa which undergoes proteolytic processing to generate 55 kDa and 42 kDa cleaved forms, and bands at the mentioned positions may be expected in Western blot application.




Images


Western Blot: PINK1 Antibody (8E10.1D6) - Azide and BSA Free [NBP2-80910] - Whole cell protein from HeLa cells treated with or without valinomycin (1 μ M, 24h) as indicated was separated by SDS-PAGE on a 7.5% polyacrylamide gel. Protein was transferred to PVDF membrane and probed with 2 μ g/ml anti-PINK1 in 1% BSA and detected with an HRP-conjugated anti-mouse secondary antibody using chemiluminescence. PINK1 is seen to be upregulated with treatment and with a molecular weight at approximately 60 kDa. Image from the standard format of this antibody.




Immunocytochemistry/Immunofluorescence: PINK1 Antibody (8E10.1D6) - Azide and BSA Free [NBP2-80910] - HeLa cells were treated with valinomycin (1 μ M, 24h) prior to being fixed in 10% buffered formalin for 10 min and permeabilized in 0.1% Triton X-100 in PBS for 10 min. Cells were incubated with NBP2-36488 at 50 μ g/ml for 1hr at room temperature, washed 3x in PBS and incubated with Alexa-DyLight-488 anti-mouse secondary antibody. PINK1 (Green) was detected at the mitochondria. Tubulin (Red) was detected using an anti-tubulin antibody with an anti-rabbit DyLight 550 secondary antibody. DNA (Blue) was counterstained with DAPI. Note: mitochondria staining might not be easily observed without treatment with valinomycin or CCCP. Image objective 40x. Image from the standard format of this antibody.



Immunohistochemistry: PINK1 Antibody (8E10.1D6) - Azide and BSA Free [NBP2-80910] - Analysis of FFPE tissue section of human hepatocellular carcinoma using PINK1 antibody (clone 8E10.1D6) at 5 μ g/ml concentration. The cancer cells developed an expected punctate to granular cytoplasmic staining with no signal in the tumor stroma. Image fr



Western Blot: PINK1 Antibody (8E10.1D6) - Azide and BSA Free [NBP2-80910] - Analysis of (A) Partial Recombinant Human PINK-1 protein with estimated molecular weight at 13kDa and (B) Human Liver lysate using PINK1 antibody clone 8E10.1D6 at 3 μ g/ml concentration, molecular weight ~64 kDa. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80910

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
BC100-494PEP	PINK1 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80910

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

