

Product Datasheet

RIPK3/RIP3 Antibody (1312D) - Azide and BSA Free NBP2-80926

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80926

Updated 2/12/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80926



NBP2-80926

RIPK3/RIP3 Antibody (1312D) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1312D
Preservative	No Preservative
Isotype	IgG
Purity	Protein A or G purified
Buffer	PBS
Target Molecular Weight	56.7 kDa

Product Description	
Host	Rabbit
Gene ID	11035
Gene Symbol	RIPK3
Species	Human
Immunogen	RIPK3/RIP3 Antibody (1312D) - Azide and BSA Free was made to a partial synthetic peptide made to an internal portion of the human RIPK3/RIP3 protein (between amino acids 470-540). This sequence has 100% homology to isoform CRA_a and isoform CRA_b.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 2 ug/ml, Immunocytochemistry/ Immunofluorescence 5 - 15 ug/ml

Images

Immunocytochemistry/Immunofluorescence: RIPK3/RIP3 Antibody (1312D) - Azide and BSA Free [NBP2-80926] - RIPK3/RIP3 was detected in immersion fixed SK-BR-3 human breast cancer cell line (positive stain) and HepG2 human human hepatocellular carcinoma cell line (negative stain) using Rabbit Monoclonal Anti-Human RIPK3/RIP3 Antibody [NBP2-76799] at 1 ug/mL for 3 hours at room temperature. Cells were stained using the NorthernLights 557-conjugated Anti-Rabbit IgG Secondary Antibody (red; [NL004]) and counterstained with DAPI (blue). Specific staining was localized to cytoplasm. Image from the standard format of this antibody.



SK-BR-3

HepG2



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80926

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00011035-P01-25ug	Recombinant Human RIPK3/RIP3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80926

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

