

Product Datasheet

Leucyl-cystinyl Aminopeptidase/LNPEP Antibody - BSA Free NBP2-81854

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81854

Updated 10/19/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81854



NBP2-81854


Leucyl-cystinyl Aminopeptidase/LNPEP Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	113 kDa
Product Description	
Host	Rabbit
Gene ID	4012
Gene Symbol	LNPEP
Species	Human
Specificity/Sensitivity	At least three isoforms are known to exist; this antibody will detect all three isoforms. Despite its predicted molecular weight, Leucyl-cystinyl Aminopeptidase/LNPEP often migrates at 165kD in SDS-PAGE.
Immunogen	LNPEP antibody was raised against a 15 amino acid synthetic peptide near the carboxy terminus of human LNPEP. The immunogen is located within amino acids 900 - 950 of LNPEP. Amino Acid Sequence: SLNGDNFRTQKLSFI
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence 20 ug/ml




Images

Western Blot: Leucyl-cystinyl Aminopeptidase/LNPEP Antibody [NBP2-81854] - Analysis of Leucyl-cystinyl Aminopeptidase/LNPEP in human lung tissue lysate with Leucyl-cystinyl Aminopeptidase/LNPEP antibody at 1 ug/ml in (A) the absence and (B) the presence of blocking peptide.



Immunocytochemistry/Immunofluorescence: Leucyl-cystinyl Aminopeptidase/LNPEP Antibody [NBP2-81854] - Immunofluorescence of Leucyl-cystinyl Aminopeptidase/LNPEP in human lung tissue with Leucyl-cystinyl Aminopeptidase/LNPEP antibody at 20 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81854

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-92071PEP	Leucyl-cystinyl Aminopeptidase/LNPEP Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81854

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

