

# Product Datasheet

## HOOK3 Antibody [PE] NBP2-81990PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-81990PE](http://www.novusbio.com/NBP2-81990PE)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-81990PE](http://www.novusbio.com/reviews/destination/NBP2-81990PE)



**NBP2-81990PE**

HOOK3 Antibody [PE]

| Product Information         |                                                                                                                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                       |
| Clonality                   | Polyclonal                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                             |
| Isotype                     | IgG                                                                                                                                                                                                            |
| Conjugate                   | PE                                                                                                                                                                                                             |
| Purity                      | Peptide affinity purified                                                                                                                                                                                      |
| Buffer                      | PBS                                                                                                                                                                                                            |
| Product Description         |                                                                                                                                                                                                                |
| Host                        | Rabbit                                                                                                                                                                                                         |
| Gene ID                     | 84376                                                                                                                                                                                                          |
| Gene Symbol                 | HOOK3                                                                                                                                                                                                          |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                              |
| Specificity/Sensitivity     | HOOK3 antibody is human, mouse and rat reactive. HOOK3 antibody is predicted to not cross-react with HOOK1 and HOOK2.                                                                                          |
| Immunogen                   | HOOK3 antibody was raised against a 19 amino acid peptide from near the amino terminus of human HOOK3. The immunogen is located within amino acids 210 - 260 of HOOK3. Amino Acid Sequence: MERLNQSDSIEDPNPAGR |
| Product Application Details |                                                                                                                                                                                                                |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry                                                                                               |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-81990PE**

---

|                    |                                              |
|--------------------|----------------------------------------------|
| NBP2-33376H        | Blue Marker Antibody (6F4-F6) [HRP]          |
| NBP2-24983         | Rabbit IgG Isotype Control [PE]              |
| H00084376-P01-10ug | Recombinant Human HOOK3 GST (N-Term) Protein |
| DKK300             | Kallikrein 3/PSA [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-81990PE](http://www.novusbio.com/reviews/submit/NBP2-81990PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

