

Product Datasheet

SUSD3 Antibody - BSA Free

NBP2-82017

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-82017

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-82017



NBP2-82017


SUSD3 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	23 kDa
Product Description	
Description	Novus Biologicals Rabbit SUSD3 Antibody - BSA Free (NBP2-82017) is a polyclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	203328
Gene Symbol	SUSD3
Species	Human, Mouse
Specificity/Sensitivity	SUSD3 antibody is human and mouse reactive. Multiple isoforms of SUSD3 are known to exist.
Immunogen	SUSD3 antibody was raised against a 16 amino acid peptide from near the center of human SUSD3. The immunogen is located within amino acids 120 - 170 of SUSD3. Amino Acid Sequence: VKKSKRRRSNRSACL
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunocytochemistry
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence 5 ug/mL, Immunocytochemistry 5 ug/ml




Images


Western Blot: SUS D3 Antibody [NBP2-82017] - Western blot analysis of SUS D3 in HeLa cell lysate with SUS D3 antibody at 1 ug/ml.




Immunocytochemistry/Immunofluorescence: SUS D3 Antibody [NBP2-82017] - Immunocytochemistry of SUS D3 in 3T3 cells with SUS D3 antibody at 5 ug/ml.



Immunocytochemistry/ Immunofluorescence: SUS D3 Antibody - BSA Free [NBP2-82017] - Immunocytochemistry of SUS D3 in 3T3 cells with SUS D3 antibody at 5 u/ml.



Immunocytochemistry/ Immunofluorescence: SUS D3 Antibody - BSA Free [NBP2-82017] - Immunocytochemistry of SUS D3 in 3T3 cells with SUS D3 antibody at 5 u/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-82017

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-82017

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

