

# Product Datasheet

## PHC3 Antibody [PerCP] NBP2-82020PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-82020PCP](http://www.novusbio.com/NBP2-82020PCP)

Updated 10/24/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-82020PCP](http://www.novusbio.com/reviews/destination/NBP2-82020PCP)



**NBP2-82020PCP**

PHC3 Antibody [PerCP]

| Product Information         |                                                                                                                                                                                                         |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                             |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                |
| Clonality                   | Polyclonal                                                                                                                                                                                              |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                      |
| Isotype                     | IgG                                                                                                                                                                                                     |
| Conjugate                   | PerCP                                                                                                                                                                                                   |
| Purity                      | Peptide affinity purified                                                                                                                                                                               |
| Buffer                      | PBS                                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                                         |
| Host                        | Rabbit                                                                                                                                                                                                  |
| Gene ID                     | 80012                                                                                                                                                                                                   |
| Gene Symbol                 | PHC3                                                                                                                                                                                                    |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                       |
| Specificity/Sensitivity     | PHC3 antibody is human, mouse and rat specific. At least six isoforms of PHC3 are known to exist.                                                                                                       |
| Immunogen                   | PHC3 antibody was raised against a 19 amino acid peptide near the amino terminus of human PHC3. The immunogen is located within amino acids 240 - 290 of PHC3. Amino Acid Sequence: KISQRDPSPESNKKGESPS |
| Product Application Details |                                                                                                                                                                                                         |
| Applications                | Western Blot, ELISA                                                                                                                                                                                     |
| Recommended Dilutions       | Western Blot, ELISA                                                                                                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-82020PCP**

---

|                    |                                                |
|--------------------|------------------------------------------------|
| NBP2-24891PCP      | Rabbit IgG Isotype Control [PerCP]             |
| H00080012-P01-10ug | Recombinant Human PHC3 GST (N-Term) Protein    |
| H00080012-T01      | PHC3 293T Cell Transient Overexpression Lysate |
| AF2400             | HNF-3 beta/FoxA2 Antibody [Unconjugated]       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-82020PCP](http://www.novusbio.com/reviews/submit/NBP2-82020PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

