

# Product Datasheet

## Human Gelsolin/GSN ELISA Kit (Colorimetric) NBP2-82209

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-82209](http://www.novusbio.com/NBP2-82209)

Updated 9/23/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-82209](http://www.novusbio.com/reviews/destination/NBP2-82209)



**NBP2-82209****Human Gelsolin/GSN ELISA Kit (Colorimetric)**


<b>Product Information</b>	
<b>Unit Size</b>	1 Kit
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Storage of components varies. See protocol for specific instructions.

<b>Product Description</b>	
<b>Gene ID</b>	2934
<b>Gene Symbol</b>	GSN
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	This kit recognizes Human Gelsolin/GSN in samples. No Significant cross-reactivity or interference between Human Gelsolin/GSN and analogues was observed.
<b>Kit Components</b>	Micro ELISA Plate, dismountable (8 wells x 12 strips); Reference Standard (2 vials); Concentrated (100X) Biotinylated Detection Ab (120 uL); Concentrated (100X) HRP Conjugate (120 uL); Reference Standard & Sample Diluent (20 mL); Biotinylated Detection Ab
<b>Standard Curve Range</b>	6.25 - 400 ng/mL
<b>Sensitivity</b>	3.75 ng/mL
<b>Inter Assay</b>	%CV < 6
<b>Intra Assay</b>	%CV < 6
<b>Assay Type</b>	Sandwich-ELISA
<b>Spiking Recovery</b>	average = 102% for serum; average = 90% for EDTA plasma; average = 95% for cell culture media
<b>Suitable Sample Type</b>	Serum, Plasma, Cell supernatant
<b>Sample Volume</b>	100 uL

<b>Product Application Details</b>	
<b>Applications</b>	ELISA
<b>Recommended Dilutions</b>	ELISA

**Images**

ELISA: Human Gelsolin/GSN ELISA Kit (Colorimetric) [NBP2-82209] - Standard Curve Reference



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-82209**

---

H00002934-P01-10ug	Recombinant Human Gelsolin/GSN GST (N-Term) Protein
202-IL-010	IL-2 [Unconjugated]
NBP1-33177	Gelsolin/GSN Antibody
AF835	Caspase-3 Antibody [Unconjugated] - Active

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-82209](http://www.novusbio.com/reviews/submit/NBP2-82209)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

