

Product Datasheet

TRPV4 Antibody - BSA Free

NBP2-82365-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-82365

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-82365



NBP2-82365-0.1ml

TRPV4 Antibody - BSA Free

| Product Information | |
|-------------------------|--|
| Unit Size | 0.1ml |
| Concentration | 0.5 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Purity | Affinity purified |
| Buffer | PBS, 2% Sucrose |
| Target Molecular Weight | 60 kDa |


| Product Description | |
|---------------------|--|
| Host | Rabbit |
| Gene ID | 59341 |
| Gene Symbol | TRPV4 |
| Species | Human |
| Immunogen | The immunogen is a synthetic peptide directed towards the N-terminal region of human TRPV4. Peptide sequence: AGPGGEVAELPGDESGTPGGGEAFPLSSLANLFEGEDGSLSPSPADASRPAG The peptide sequence for this immunogen was taken from within the described region. |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin |



Images


Western Blot: TRPV4 Antibody [NBP2-82365] - Host: Rabbit. Target Name: TRPV4. Sample Type: HT1080 Whole cell lysates. Antibody Dilution: 1.0ug/ml



Immunohistochemistry-Paraffin: TRPV4 Antibody [NBP2-82365] - Rabbit Anti-TRPV4 Antibody. Formalin Fixed Paraffin Embedded Tissue: Human Adult heart. Observed Staining: Cytoplasmic. Primary Antibody Concentration: 1:100. Secondary Antibody: Donkey anti-Rabbit-Cy2/3. Secondary Antibody Concentration: 1:200. Magnifica



Western Blot: TRPV4 Antibody [NBP2-82365] - WB Suggested Anti-TRPV4 antibody Titration: 1 ug/mL. Sample Type: Human heart





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-82365-0.1ml

| | |
|-------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBL1-17343 | TRPV4 Overexpression Lysate |
| NB110-40763 | TRPA1 Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-82365

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

