

Product Datasheet

LEPROTL1 Antibody - BSA Free

NBP2-83143-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-83143

Updated 2/26/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-83143



NBP2-83143-0.1ml

LEPROTL1 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Product Description	
Host	Rabbit
Gene ID	23484
Gene Symbol	LEPROTL1
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human LEPROTL1. Peptide sequence: FVLFFYILSPIPYCIARRLVDDTDAMSNACKELAIFLTGGIVVSAFGLPI The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



Images

Western Blot: LEPROTL1 Antibody [NBP2-83143] - Host: Rabbit. Target Name: LEPROTL1. Sample Type: Human HepG2. Antibody Dilution: 1.0ug/ml LEPROTL1 is strongly supported by BioGPS gene expression data to be expressed in Human HepG2 cells

LEPROTL1




Western Blot: LEPROTL1 Antibody [NBP2-83143] - Host: Rabbit. Target Name: LEPROTL1. Sample Type: Fetal Liver lysates. Antibody Dilution: 1.0ug/ml




Western Blot: LEPROTL1 Antibody [NBP2-83143] - Host: Rabbit. Target Name: LEPROTL1. Sample Type: Human HeLa. Antibody Dilution: 1.0ug/ml LEPROTL1 is supported by BioGPS gene expression data to be expressed in HeLa

LEPROTL1



Western Blot: LEPROTL1 Antibody [NBP2-83143] - Host: Rabbit. Target Name: LEPROTL1. Sample Type: Human 721_B. Antibody Dilution: 1.0ug/ml LEPROTL1 is supported by BioGPS gene expression data to be expressed in 721_B

LEPROTL1





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-83143-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-14626PEP	LEPROTL1 Recombinant Protein Antigen
291-G1-200	IGF-I/IGF-1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-83143

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

