

Product Datasheet

RPS28 Antibody - BSA Free

NBP2-85673-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-85673

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-85673



NBP2-85673-0.1ml

RPS28 Antibody - BSA Free

Product Information	
Unit Size	0.1ml
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Product Description	
Description	Novus Biologicals Rabbit RPS28 Antibody - BSA Free (NBP2-85673) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6234
Gene Symbol	RPS28
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human RPS28. Peptide sequence: RTGSQGQCTQVRVEFMDDTSRSIIRNVKGPVREGDVLTLLESEREARRLR The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin



Images


Western Blot: RPS28 Antibody [NBP2-85673] - WB Suggested Anti-RPS28 Antibody. Titration: 1.0 ug/ml. Positive Control: MDA-MB-435S Whole Cell RPS28 is strongly supported by BioGPS gene expression data to be expressed in Human MDA-MB435 cells



Immunohistochemistry-Paraffin: RPS28 Antibody [NBP2-85673] - Rabbit Anti-RPS28 Antibody. Formalin Fixed Paraffin Embedded Tissue: Human Adult Liver. Observed Staining: Cytoplasm in Kupffer cells only, strong signal, very low tissue distribution. Primary Antibody Concentration: 1:100. Secondary Antibody: Donkey anti-Rabbit-Cy3. Secondary Antibody Concentration: 1:200. Magnification: 20X.



Western Blot: RPS28 Antibody [NBP2-85673] - Host: Rabbit. Target: RPS28. Positive control (+): THP-1 (N30). Negative control (-): A549 (N03). Antibody concentration: 3ug/ml



Host: Rabbit
Target: name: RPS28
Positive control: ~25ug THP1 Cell lysate (THP1)
Negative control: ~25ug A549 Cell lysate (A549)
Antibody concentration: 3ug/ml

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-85673-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP1-92353PEP	RPS28 Recombinant Protein Antigen
NB500-171	Histone H3 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-85673

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

