

Product Datasheet

CLK4 Antibody - BSA Free

NBP2-87195

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-87195

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-87195



NBP2-87195

CLK4 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Product Description	
Description	Novus Biologicals Rabbit CLK4 Antibody - BSA Free (NBP2-87195) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	57396
Gene Symbol	CLK4
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of human CLK4. Peptide sequence: RSRRSSPKRKRNRHCSSHQSRSKSHRRKRSR SIEDDEEGHLICQSGDVLR The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



Images

Western Blot: CLK4 Antibody [NBP2-87195] - Sample type: human microvascular endothelial cells (25ug). Primary Dilution: 1:1000. Secondary Antibody: Goat anti-Rabbit-HRP. Secondary Dilution: 1:5000. Image Submitted by: Andreas Eisenreich. Charite Universitatsmedizin Berlin

CLK4




See Immunoblot 2 Data for more information.

Western Blot: CLK4 Antibody [NBP2-87195] - Host: Rabbit. Target Name: CLK4. Sample Tissue: Human PANC1 Whole Cell. Antibody Dilution: 1ug/ml



Host: Rabbit
Target Name: CLK4
Sample Type: PANC1 Cell Lysate
Antibody Dilution: 1.0µg/ml

Western Blot: CLK4 Antibody [NBP2-87195] - WB Suggested Anti-CLK4 Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:62500. Positive Control: Human brain



Western Blot: CLK4 Antibody [NBP2-87195] - Host: Rabbit. Target: CLK4. Positive control (+): Human lung (LU). Negative control (-): HEK293 (HEK293). Antibody concentration: 1ug/ml



Host: Rabbit
Target name: CLK4
Positive control: ~25µg Human Lung [LU]
Negative control: ~25µg HEK293 Cell Lysate (2T)
Antibody concentration: 1ug/ml



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-87195

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-47505PEP	CLK4 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-87195

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

