

Product Datasheet

DPF3 Antibody - BSA Free NBP2-87297

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-87297

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-87297



NBP2-87297

DPF3 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose

Product Description


Description	Novus Biologicals Rabbit DPF3 Antibody - BSA Free (NBP2-87297) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8110
Gene Symbol	DPF3
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human DPF3. Peptide sequence: KDGFTSESTTLEALLRGEGVEKKVDAREEESIQEIQRVLENDENVEEGNE The peptide sequence for this immunogen was taken from within the described region.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml

Images


Western Blot: DPF3 Antibody [NBP2-87297] - Host: Rabbit. Target Name: DPF3. Sample Type: NCI-H226. Antibody Dilution: 1.0ug/ml DPF3 is supported by BioGPS gene expression data to be expressed in NCIH226

DPF3


Western Blot: DPF3 Antibody [NBP2-87297] - WB Suggested Anti-DPF3 Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:12500. Positive Control: MCF7 cell lysate




Western Blot: DPF3 Antibody [NBP2-87297] - Host: Rabbit. Target Name: DPF3. Sample Type: MCF7. Antibody Dilution: 1.0ug/ml



Western Blot: DPF3 Antibody [NBP2-87297] - Host: Rabbit. Target Name: DPF3. Sample Type: HT1080. Antibody Dilution: 1.0ug/ml



Western Blot: DPF3 Antibody [NBP2-87297] - Host: Rabbit. Target Name: DPF3. Sample Tissue: Human HepG2 Whole Cell. Antibody Dilution: 1ug/ml





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-87297

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-57317PEP	DPF3 Recombinant Protein Antigen
NBP2-29661-1Kit	Human, Mouse, Rat RelA/NFkB p65 ELISA Kit (Colorimetric)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-87297

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

