

Product Datasheet

POU2F3 Antibody - BSA Free

NBP2-88075

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-88075

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-88075



NBP2-88075


POU2F3 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Protein A purified
Buffer	PBS, 2% Sucrose
Product Description	
Description	Novus Biologicals Rabbit POU2F3 Antibody - BSA Free (NBP2-88075) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	25833
Gene Symbol	POU2F3
Species	Human, Mouse, Rat
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human POU2F3. Peptide sequence: KMSGDVADSTDARSTLSQVEPGNDRKGLDFNRQIKTEDLSDSLQQTLSHR The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin




Images

Western Blot: POU2F3 Antibody [NBP2-88075] - Host: Rat. Target Name: POU2F3. Sample Tissue: Rat Brain. Antibody Dilution: 1ug/ml



Immunohistochemistry: POU2F3 Antibody [NBP2-88075] - Sample Type: Mouse tongue tissue. Primary Antibody Dilution: 1:100. Secondary Antibody: Anti-rabbit-Cy3. Secondary Antibody Dilution: 1:500. Color/Signal Descriptions: Red: POU2F3. Gene Name: POU2F3. Submitted by: Dr. Hong Wang, Monell Chemical Senses Cen



See IHC 1 Data and Customer Feedback for more information

Western Blot: POU2F3 Antibody [NBP2-88075] - WB Suggested Anti-POU2F3 Antibody. Titration: 1.25 ug/ml. Positive Control: HepG2 Whole Cell





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-88075

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP1-83966PEP	POU2F3 Recombinant Protein Antigen
NBL1-14618	POU2F3 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-88075

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

