

Product Datasheet

Human Tmprss2 ELISA Kit (Colorimetric) NBP2-89170

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-89170

Updated 1/8/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-89170




NBP2-89170

Human TMPRSS2 ELISA Kit (Colorimetric)


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	7113
Gene Symbol	TMPRSS2
Species	Human
Specificity/Sensitivity	This Human TMPRSS2 ELISA Kit (Colorimetric) recognizes Human TMPRSS2 in samples. No Significant cross-reactivity or interference between Human TMPRSS2 and analogues was observed.
Kit Components	Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Manual
Standard Curve Range	0.16-10 ng/mL
Sensitivity	0.09 ng/mL
Inter Assay	CV% < 5.0%
Intra Assay	CV% < 4.5%
Assay Type	Sandwich-ELISA
Spiking Recovery	89-108%
Suitable Sample Type	Serum, plasma, cell supernatant, and other biological fluids
Sample Volume	100 ul
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human Tmprss2 ELISA Kit (Colorimetric) [NBP2-89170] - Samples were spiked with high concentrations of Human Tmprss2 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human Tmprss2 ELISA Kit (Colorimetric) [NBP2-89170] - Standard Curve Reference



Publications

Tampe D, Hakroush S, BOsherz M et al. Urinary Levels of SARS-CoV-2 Nucleocapsid Protein Associate With Risk of AKI and COVID-19 Severity: A Single-Center Observational Study *Frontiers in Medicine* 2021-05-25 [PMID: 34113632] (ELISA)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-89170

H00007113-Q01-10ug	Recombinant Human TMPRSS2 GST (N-Term) Protein
AF933	ACE-2 Antibody [Unconjugated]
NBP1-20984	TMPRSS2 Antibody
AF262	Amphiregulin Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-89170

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

