

Product Datasheet

Serpin A8/Angiotensinogen Antibody (1) - Azide and BSA Free NBP2-89736-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-89736

Updated 7/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-89736



NBP2-89736-100ul

Serpine A8/Angiotensinogen Antibody (1) - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 um filtered solution in PBS


Product Description	
Description	This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Gene ID	183
Gene Symbol	AGT
Species	Human
Specificity/Sensitivity	No cross-reactivity in ELISA with: Human SerpinA1 / A1AT Human SerpinA3 / SERPINA3 / AACT Human SerpinA5 / Protein C inhibitor Human SerpinD1 / HC2 Human SerpinF2 / SERPINF2 Human SerpinG1 / C1IN Human SerpinI1 / Protease inhibitor 12 Human SerpinC1 /SERPINC1/ Antithrombin-III
Immunogen	This antibody was obtained from a rabbit immunized with purified, recombinant Human Serpin A8/Angiotensinogen (Accession#: NP_000020.1; Met 1-Ala 485).

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:2000, Flow Cytometry 1:25-1:100, ELISA 1:5000-1:10000, Immunocytochemistry/ Immunofluorescence 1:20-1:100, Immunoprecipitation 4-6 uL/mg of lysate




Images


Western Blot: Serpin A8/Angiotensinogen Antibody (1) [NBP2-89736] -
 Lane A: HepG2 Whole Cell Lysate Lysates/proteins at 30 ug per lane.
 Secondary Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.
 Developed using the Odyssey technique. Performed under reducing
 conditions. Predicted band size:53 kDa Observed band size:53 kDa




Flow Cytometry: Serpin A8/Angiotensinogen Antibody (1) [NBP2-89736]
 - Analysis of Human AGT expression on HepG2 cells. The cells were
 treated ,stained with purified anti-Human AGT, then a FITC-conjugated
 second step antibody. The fluorescence histograms were derived from
 gated events with the forward and side light-scatter characteristics of
 intact cells.



Immunoprecipitation: Serpin A8/Angiotensinogen Antibody (1) [NBP2-89736] - Lane A:0.5 mg HepG2 Whole Cell Lysate 2 uL anti-SerpinA8 rabbit monoclonal antibody and 15 ul of 50 % Protein G agarose.
 Primary antibody: Anti-SerpinA8 rabbit monoclonal antibody,at 1:200 dilution
 Secondary antibody: Clean-Bloto IP Detection Reagent (HRP) at 1:500 dilution
 Developed using the DAB staining technique. Performed under reducing conditions. Predicted band size: 51 kDa Observed band size: 51 kDa



Immunocytochemistry/Immunofluorescence: Serpin A8/Angiotensinogen Antibody (1) [NBP2-89736] - Staining of Human AGT in Hela or SKBR3 cells. Cells were fixed with 4% PFA, permeabilized with 1% Triton X-100 in PBS, blocked with 10% serum, and incubated with NBP2-89736 antibody (1:60). Then cells were stained with the Alexa Fluor(R) 488-conjugated Goat Anti-rabbit IgG secondary antibody (left panel, captured by laser confocal scanning microscope; right panel, captured by fluorescence microscope), countstained with DAPI (blue). Positive staining was localized to cytoplasm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-89736-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-30027PEP	Serpin A8/Angiotensinogen Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-89736

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

