

Product Datasheet

DGK-theta Antibody - BSA Free NBP2-92780-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-92780

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-92780



NBP2-92780-0.1ml

DGK-theta Antibody - BSA Free


| Product Information | |
|-------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | PBS with 50% glycerol, pH7.3. |
| Target Molecular Weight | 101 kDa |

| Product Description | |
|---------------------|--|
| Host | Rabbit |
| Gene ID | 1609 |
| Gene Symbol | DGKQ |
| Species | Human, Mouse |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 683-942 of human DGKQ (NP_001338.2). GRVLRWGAGYSGEDPFSVLLSVDEADAVLMDRWTILLDAHEAGSAENDTADA EPPKIVQMSNYCGIGIDAELSLDFHQAREEPEGKFTSRLHNKGVYVRVGLQKIS HSRSLHKQIRLQVERQEVELPSIEGLIFINIPSWGSGADLWGSDDSTRFEKPRM DDGLLEVVGVTGVVHMGQVQGGLRSGIRIAQGSYFRVTLKATPVQVDGEPW VQAPGHMIIISAAGPKVHMLRKAKQKPRRAGTTRDARADAAPAPESDPR |

| Product Application Details | |
|-----------------------------|-----------------------------|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot 1:200 - 1:2000 |

Images

Western Blot: DGK-theta Antibody [NBP2-92780] - Analysis of extracts of various cell lines, using DGK-theta at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-92780-0.1ml

| | |
|---------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP1-85317PEP | DGK-theta Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-92780

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

