

Product Datasheet

Cochlin Antibody - BSA Free

NBP2-92904-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-92904

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-92904



NBP2-92904-0.1ml


Cochlin Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	59 kDa
Product Description	
Description	Novus Biologicals Rabbit Cochlin Antibody - BSA Free (NBP2-92904) is a polyclonal antibody validated for use in IHC, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1690
Gene Symbol	COCH
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human COCH (NP_001128530.1). GPAGSEGAAPAIATCFTRGLDIRKEKADVLCPPGGCPLEEFVYGNIVYASVSSIC GAAVHRGVISNSGGPVRVYSLPGRENYSSVDANGIQSQMLSRWSASFVTGK KSSTQEATGQAVSTAHPPTGKRLKKTPEKKTGNKDCADIAFLIDGSFNIGQRR FNLQKNFVGKVALMLGIGTEGPHVGLVQASEHPKIEFYLNFTSAKDVLFAIKEV GFRGGNSNTGKALKHTAQKFFTVDA
Product Application Details	
Applications	Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	ELISA, Immunohistochemistry 1:50-1:100, Immunocytochemistry/ Immunofluorescence 1:50-1:100, Immunohistochemistry-Paraffin




Images

Immunohistochemistry-Paraffin: Cochlin Antibody [NBP2-92904] - Paraffin-embedded human prostate using Cochlin .



Immunocytochemistry/ Immunofluorescence: Cochlin Antibody - BSA Free [NBP2-92904] - Immunofluorescence analysis of HeLa cells using Cochlin Rabbit pAb . Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H +L) at 1:500 dilution. Blue: DAPI for nuclear staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-92904-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP3-18160	Recombinant Mouse Cochlin FLAG Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-92904

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

