

# Product Datasheet

## **CARD11/CARMA1 Antibody - BSA Free** **NBP2-92910-0.1ml**

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-92910](http://www.novusbio.com/NBP2-92910)

Updated 11/7/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-92910](http://www.novusbio.com/reviews/destination/NBP2-92910)



**NBP2-92910-0.1ml**


CARD11/CARMA1 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                  |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                      |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                       |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                              |
| Purity                      | Affinity purified                                                                                                                                                                                                                                                                                                |
| Buffer                      | PBS with 50% glycerol, pH7.3.                                                                                                                                                                                                                                                                                    |
| Target Molecular Weight     | 133 kDa                                                                                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                                                                                  |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 84433                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | CARD11                                                                                                                                                                                                                                                                                                           |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                |
| Immunogen                   | Recombinant fusion protein containing a sequence corresponding to amino acids 985-1154 of human CARD11 (NP_115791.3).<br>AKTLVQRLLNSGGAMEFTICKSDIVTRDEFLRRQKTETIIYSREKNPNAFECIAPA<br>NIEAVALAKNKHCLLEAGIGCTRDLIKSNIYPIVLFIRVCEKNIKRFRKLLPRPETEE<br>EFLRVCRLKEKELEALPCLYATVEPDMWGSVEELLRVVKDKIGEEQRKTIWVD<br>EDQL |
| Product Application Details |                                                                                                                                                                                                                                                                                                                  |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                                                                                                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot 1:500 - 1:2000, Immunocytochemistry/ Immunofluorescence 1:50 - 1:100                                                                                                                                                                                                                                |




## Images

Western Blot: CARD11/CARMA1 Antibody [NBP2-92910] - Analysis of extracts of various cell lines, using CARD11/CARMA1 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.



Immunocytochemistry/Immunofluorescence: CARD11/CARMA1 Antibody [NBP2-92910] - Immunofluorescence analysis of A549 cells using CARD11/CARMA1 antibody (NBP2-92910). Blue: DAPI for nuclear staining.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-92910-0.1ml**

---

|                    |                                                      |
|--------------------|------------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]        |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]  |
| NBP2-24891         | Rabbit IgG Isotype Control                           |
| H00084433-Q01-10ug | Recombinant Human CARD11/CARMA1 GST (N-Term) Protein |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-92910](http://www.novusbio.com/reviews/submit/NBP2-92910)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

