

# Product Datasheet

## CSRP2BP Antibody - BSA Free

### NBP2-92912-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-92912](http://www.novusbio.com/NBP2-92912)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-92912](http://www.novusbio.com/reviews/destination/NBP2-92912)



**NBP2-92912-0.1ml**

CSRP2BP Antibody - BSA Free

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| Clonality           | Polyclonal                                                                                  |
| Preservative        | 0.02% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Purity              | Affinity purified                                                                           |
| Buffer              | PBS (pH 7.3), 50% glycerol                                                                  |


| Product Description |                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit CSRP2BP Antibody - BSA Free (NBP2-92912) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                         |
| Host                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                     |
| Gene ID             | 57325                                                                                                                                                                                                                                                                                                                                                                                                      |
| Gene Symbol         | KAT14                                                                                                                                                                                                                                                                                                                                                                                                      |
| Species             | Human, Rat                                                                                                                                                                                                                                                                                                                                                                                                 |
| Immunogen           | Recombinant fusion protein containing a sequence corresponding to amino acids 523-782 of human CSRP2BP (NP_065397.3).<br>DGIYGAKEGGISRLPAGQATYRTTCQDFRILDYQTSLSRKGFRHQTTKFLYR<br>LVGSEDMAVDQSIVSPYTSRILKPYIRRDYETKPPKLQLLSQIRSHLHRSDPHWT<br>PEPDAPLDYCYVRPNHIPTINSMCQEFFWPGIDLSECLQYPDFSVVLYKKVIA<br>FGFMVPDVKYNEAYISFLFVHPEWRRAGIATFMIYHLIQTCMGKDVTLHVSASN<br>PAMLLYQKFGFKTEEYVLDFYDKYYPLESTECKHAFFLRRLRR |

| Product Application Details |                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                              |
| Recommended Dilutions       | Western Blot 1:1000 - 1:5000, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200 |



## Images

Western Blot: CSRP2BP Antibody [NBP2-92912] - Analysis of various lysates, using CSRP2BP Rabbit pAb at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunocytochemistry/ Immunofluorescence: CSRP2BP Antibody - BSA Free [NBP2-92912] - Immunofluorescence analysis of A549 cells using CSRP2BP antibody . Blue: DAPI for nuclear staining.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-92912-0.1ml**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-57876PEP | CSRP2BP Recombinant Protein Antigen                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-92912](http://www.novusbio.com/reviews/submit/NBP2-92912)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

