

Product Datasheet

Proteasome 20S alpha 3 Antibody - BSA Free NBP2-94032-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-94032

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-94032



NBP2-94032-0.1ml

Proteasome 20S alpha 3 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol

Product Description	
Description	Novus Biologicals Rabbit Proteasome 20S alpha 3 Antibody - BSA Free (NBP2-94032) is a polyclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5684
Gene Symbol	PSMA3
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human PSMA3 (NP_002779.1). MSSIGTGYDLSASTFSPDGRVFQVEYAMKAVENTSSTAIGIRCKDGVVFGVEKLV LSKLYEEGSNKRLFNVDHRHVGMAVAGLLADARSLADIAREEASNFRSNGFYINIP LKHLADRVAMYVHAYTLYSVRPFGCSFMLGSYSVNDGAQLYMIDPSGVSYG YWGCAIGKARQAAKTEIEKLQMKEMTCRDIVKEVAKIIVHDEVKDKAFELELS WVGELTNGRHEIVPKDIREEAKEYAKESLKEEDESDDDNM


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunoprecipitation 1:50 - 1:100, Immunohistochemistry-Paraffin

Images


Western Blot: Proteasome 20S alpha 3 Antibody [NBP2-94032] - Analysis of extracts of various cell lines, using PSMA3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 3 Antibody [NBP2-94032] - Rat brain using PSMA3 Rabbit pAb at dilution of 1:100 (40x lens).




Western Blot: Proteasome 20S alpha 3 Antibody [NBP2-94032] - Analysis of 200ug extracts of HL-60 cells, using Proteasome 20S alpha 3 at a dilution of 1:1000.




Immunohistochemistry-Paraffin: Proteasome 20S alpha 3 Antibody [NBP2-94032] - Human colon carcinoma using PSMA3 Rabbit pAb at dilution of 1:100 (40x lens).



Immunohistochemistry-Paraffin: Proteasome 20S alpha 3 Antibody [NBP2-94032] - Mouse heart using PSMA3 Rabbit pAb at dilution of 1:100 (40x lens).



Immunocytochemistry/ Immunofluorescence: Proteasome 20S alpha 3
Antibody - Azide and BSA Free [Proteasome 20S alpha 3] -
Immunofluorescence analysis of Mouse kidney tissue using Proteasome
20S alpha 3 Rabbit pAb at a dilution of 1:400 (40x lens). Secondary
antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) at 1:500 dilution.
Blue: DAPI for nuclear staining. High pressure antigen retrieval
performed with 0.01M Citrate Buffer (pH 6.0) prior to IF staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-94032-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-72526-100ug	Recombinant Human Proteasome 20S alpha 3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-94032

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

