

Product Datasheet

RPL36 Antibody - BSA Free

NBP2-94161-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-94161

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-94161



NBP2-94161-0.1ml

RPL36 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	12 kDa

Product Description	
Host	Rabbit
Gene ID	25873
Gene Symbol	RPL36
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human RPL36 (NP_056229.2). MALRYPMAGVLNKGHKVTKNVSKPRHSRRRGRLTKHTKFVRDMIREVCGFAP YERRAMELLKVSCKDKRALKFIIKRVGTHIRAKRKREELSNVLAAMRKAANKD


Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry 1:50-1:200, Immunocytochemistry/Immunofluorescence 1:50 - 1:200, Immunohistochemistry-Paraffin

Images


Western Blot: RPL36 Antibody [NBP2-94161] - Analysis of extracts of various cell lines, using RPL36 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit . Exposure time: 10s.




Immunocytochemistry/Immunofluorescence: RPL36 Antibody [NBP2-94161] - Confocal immunofluorescence analysis of HeLa cells using RPL36 Polyclonal Antibody (NBP2-94161) at dilution of 1:200. Blue: DAPI for nuclear staining.




Immunohistochemistry-Paraffin: RPL36 Antibody [NBP2-94161] - Immunohistochemistry of paraffin-embedded mouse heart using RPL36 Rabbit pAb (NBP2-94161) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunocytochemistry/Immunofluorescence: RPL36 Antibody [NBP2-94161] - Analysis of U2OS cells using RPL36 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunocytochemistry/Immunofluorescence: RPL36 Antibody [NBP2-94161] - Confocal immunofluorescence analysis of HeLa cells using RPL36 antibody (NBP2-94161) at dilution of 1:200. Blue: DAPI for nuclear staining.




Immunohistochemistry-Paraffin: RPL36 Antibody [NBP2-94161] - Paraffin-embedded mouse brain using RPL36 .



Immunohistochemistry-Paraffin: RPL36 Antibody [NBP2-94161] - Immunohistochemistry of paraffin-embedded rat brain using RPL36 Rabbit pAb (NBP2-94161) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry-Paraffin: RPL36 Antibody [NBP2-94161] - Immunohistochemistry of paraffin-embedded human liver cancer using RPL36 Rabbit pAb (NBP2-94161) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-94161-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00025873-P01-10ug	Recombinant Human RPL36 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-94161

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

