

# Product Datasheet

## PFDN4 Antibody [DyLight 405]

### NBP3-00251V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-00251V](http://www.novusbio.com/NBP3-00251V)

Updated 10/17/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-00251V](http://www.novusbio.com/reviews/destination/NBP3-00251V)



**NBP3-00251V**

PFDN4 Antibody [DyLight 405]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.          |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Polyclonal   |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG  |
| Conjugate                   | DyLight 405  |
| Purity                      | Antigen and protein A Affinity-purified  |
| Buffer                      | 50mM Sodium Borate   |
| Product Description         |  |
| Host                        | Rabbit   |
| Gene ID                     | 5203   |
| Gene Symbol                 | PFDN4  |
| Species                     | Human  |
| Immunogen                   | Produced in rabbits immunized with purified, recombinant Human PFDN4 (Uniprot#: Q9NQP4; Met1-Ser134) |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                    |
| Product Application Details |  |
| Applications                | ELISA, Immunohistochemistry-Paraffin   |
| Recommended Dilutions       | ELISA, Immunohistochemistry-Paraffin   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-00251V**

---

|                   |  |
|-------------------|--|
| NBP2-24891V       | Rabbit IgG Isotype Control [DyLight 405] |
| NBP1-50973-0.05mg | Recombinant Human PFDN4 His Protein      |
| DY2037            | Complement Component C5a [Biotin]        |
| NBL1-14308        | PFDN4 Overexpression Lysate              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-00251V](http://www.novusbio.com/reviews/submit/NBP3-00251V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

