

Product Datasheet

Human Phospholipid Scramblase 1/PLSCR1 ELISA Kit (Colorimetric) NBP3-00413

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-00413

Updated 1/8/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-00413




NBP3-00413**Human Phospholipid Scramblase 1/PLSCR1 ELISA Kit (Colorimetric)**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	5359
Gene Symbol	PLSCR1
Species	Human
Kit Components	Reference Standard, Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Substrate Reagent, Stop Solution, Plate Sealer, Concentrated Biotinylated Detection Ab (100X), Concentrated HRP Conjugate (100X), Concentrated Wash Buffer (25X), Micro ELISA Plate (Dismountable)
Standard Curve Range	0.31 - 20 ng/mL
Sensitivity	0.19 ng/mL
Inter Assay	%CV < 5.77
Intra Assay	%CV < 5.58
Assay Type	Sandwich-ELISA
Spiking Recovery	90 - 103%
Suitable Sample Type	Serum, plasma, cell supernatant, and other biological fluids
Sample Volume	100 uL
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Human Phospholipid Scramblase 1/PLSCR1 ELISA Kit (Colorimetric) [NBP3-00413] - Samples were spiked with high concentrations of Human Phospholipid Scramblase 1/PLSCR1 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human Phospholipid Scramblase 1/PLSCR1 ELISA Kit (Colorimetric) [NBP3-00413] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-00413

NBP2-23329	Recombinant Human Phospholipid Scramblase 1/PLSCR1 His Protein
236-EG-200	EGF [Unconjugated]
NBP1-32588	Phospholipid Scramblase 1/PLSCR1 Antibody
AF835	Caspase-3 Antibody [Unconjugated] - Active

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-00413

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

