

Product Datasheet

Human WISP-1/CCN4 ELISA Kit (Colorimetric) NBP3-00420

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-00420

Updated 1/8/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-00420




NBP3-00420**Human WISP-1/CCN4 ELISA Kit (Colorimetric)**


| Product Information | |
|------------------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Storage of components varies. See protocol for specific instructions. |
| Product Description | |
| Gene ID | 8840 |
| Gene Symbol | CCN4 |
| Species | Human |
| Specificity/Sensitivity | This kit recognizes Human WISP-1/CCN4 in samples. No significant cross-reactivity or interference between Human WISP-1/CCN4 and analogues was observed. |
| Kit Components | Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Manual |
| Standard Curve Range | 78.13 - 5000 pg/mL |
| Sensitivity | 46.88 pg/mL |
| Inter Assay | 5.46% |
| Intra Assay | 5.12% |
| Assay Type | Sandwich-ELISA |
| Spiking Recovery | 87 - 106% |
| Suitable Sample Type | Serum, plasma, cell supernatant, and other biological fluids |
| Sample Volume | 100 ul |
| Product Application Details | |
| Applications | ELISA |
| Recommended Dilutions | ELISA |

Images

ELISA: Human WISP-1/CCN4 ELISA Kit (Colorimetric) [NBP3-00420] - Samples were spiked with high concentrations of Human WISP-1/CCN4 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human WISP-1/CCN4 ELISA Kit (Colorimetric) [NBP3-00420] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-00420

| | |
|---------------|---|
| NBP1-86570PEP | WISP-1/CCN4 Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| AF1680 | WISP-1/CCN4 Antibody |
| 355-BM-010 | BMP-2 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-00420

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

