

# Product Datasheet

## Frequenin Antibody - BSA Free NBP3-03315-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-03315](http://www.novusbio.com/NBP3-03315)

Updated 11/7/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-03315](http://www.novusbio.com/reviews/destination/NBP3-03315)



**NBP3-03315-100ul**

Frequenin Antibody - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | PBS with 50% glycerol, pH7.3.                                                               |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                 |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                          |
| <b>Gene ID</b>     | 23413                                                                                                                                                                                                                                                                                                                           |
| <b>Gene Symbol</b> | NCS1                                                                                                                                                                                                                                                                                                                            |
| <b>Species</b>     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen</b>   | Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human Frequenin (NP_055101.2).<br>MGKSNSKLKPEVVEELTRKTYFTEKEVQQWYKGFIDCPGQLDAAGFQKIYK<br>QFFPFGDPTKFATFVFNVDENKDGRIEFSEFIQALSVTSRGTLDEKLRWAFKL<br>YDLNDGYITRNEMLDIVDAIYQMVGNVELPEEENTPEKRVDRIFAMMDKNAD<br>GKLTQLQEFQEGSKADPSIVQALSLYDGLV |

**Product Application Details**

|                              |                           |
|------------------------------|---------------------------|
| <b>Applications</b>          | Western Blot              |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:2000 |

**Images**

Immunohistochemistry-Paraffin: Frequenin Antibody [NBP3-03315] - Analysis of extracts of various cell lines, using Frequenin antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL B





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-03315-100ul**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control                          |
| NBP1-44484 | Recombinant Human Frequentin His Protein            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-03315](http://www.novusbio.com/reviews/submit/NBP3-03315)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

