

Product Datasheet

AXUD1 Antibody - BSA Free

NBP3-03327-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-03327

Updated 4/20/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-03327



NBP3-03327-100ul


AXUD1 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	64 kDa
Product Description	
Host	Rabbit
Gene ID	64651
Gene Symbol	CSRNP1
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human AXUD1 (NP_149016.2). MTGLLKRKFDQLDEDNSSVSSSSSSSGCQSRSCSPSSSVSRAWDSEEEGPW DQMPLPDRDFCGPRSFTPLSILKRARRERPGRVAFDGITVFYFPRCQGFTSVP SRGGCTLGMALRHSACRRFSLAEFAQEQARARHEKLRQRLKEEKLEMLQWKL SAAGVPQAEAGLPPVDAIDDASVEEDLAVAVAGGRLEEVSFLQ
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:2000, ELISA, Immunohistochemistry 1:50-1:200, Immunocytochemistry/ Immunofluorescence 1:50-1:100, Immunohistochemistry-Paraffin




Images


Immunocytochemistry/Immunofluorescence: AXUD1 Antibody [NBP3-03327] - Analysis of A549 cells using AxUD1 antibody .



Immunohistochemistry-Paraffin: AXUD1 Antibody [NBP3-03327] - Paraffin-embedded rat kidney using AxUD1 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry-Paraffin: AXUD1 Antibody [NBP3-03327] - Paraffin-embedded mouse brain using AxUD1 antibody at dilution of 1:200 (40x lens).




Western Blot: AXUD1 Antibody - BSA Free [AXUD1] - Western blot analysis of various lysates, using AXUD1 Rabbit pAb at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit .

Exposure time: 10s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-03327-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00064651-P01-10ug	Recombinant Human AXUD1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-03327

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

