

Product Datasheet

ARL4 Antibody - Azide and BSA Free NBP3-03356-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-03356

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-03356



NBP3-03356-100ul

ARL4 Antibody - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	22 kDa

Product Description


Description	Novus Biologicals Rabbit ARL4 Antibody - Azide and BSA Free (NBP3-03356) is a polyclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10124
Gene Symbol	ARL4A
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-200 of human ARL4 (NP_001032241.1). ISENQGVPLIVANKQDLRNLSLSEIEKLLAMGELSSSTPWHLQPTCAIIGDGL KEGLEKLHDMIKRRKMLRQQKKKR

Product Application Details


Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:50-1:200, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin

Images


Immunohistochemistry-Paraffin: ARL4 Antibody [NBP3-03356] - Paraffin-embedded human liver cancer using ARL4 antibody at dilution of 1:100 (40x lens).




Immunocytochemistry/Immunofluorescence: ARL4 Antibody [NBP3-03356] - Analysis of NIH-3T3 cells using ARL4 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry-Paraffin: ARL4 Antibody [NBP3-03356] - Paraffin-embedded rat liver using ARL4 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry-Paraffin: ARL4 Antibody [NBP3-03356] - Paraffin-embedded mouse heart using ARL4 antibody at dilution of 1:100 (40x lens).



Immunocytochemistry/ Immunofluorescence: ARL4 Antibody - Azide and BSA Free [ARL4] - Immunofluorescence analysis of U-2 OS cells using ARL4 Rabbit pAb at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-03356-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-72396-50ug	Recombinant Human ARL4 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-03356

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

