

Product Datasheet

GABA-A R delta Antibody - Azide and BSA Free NBP3-03679-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-03679

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-03679



NBP3-03679-100ul

GABA-A R delta Antibody - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Proclin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	51 kDa

Product Description

Description	Novus Biologicals Rabbit GABA-A R delta Antibody - Azide and BSA Free (NBP3-03679) is a polyclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2563
Gene Symbol	GABRD
Species	Human, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-248 of human GABA-A R delta (NP_000806.2)). YVGSNLEISWLPNLDGLIAGYARNFRPGIGGPPVNVALALEVASIDHISEANMEY TMTVFLHQSWRDSRLSYNHTNETLGLDSRFVDKLWLPDTFIVNAKSAWFHDVT VENKLIRLQPDGVILYSIRITSTVACDMDLAKYPMDEQECMLDLESYGYSSEDIV YYWSESQEHIHGLDKLQLAQFTITSYRFTTELMNFKSAGQFPRLSLHFHLRRNR G

Product Application Details

Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500-1:2000, ELISA Recommended starting concentration is 1 µg/mL

Images

Western Blot: GABA-A R delta Antibody - Azide and BSA Free [GABA-A R delta] - Western blot analysis of lysates from MCF7 cells, using GABA-A R delta Rabbit pAb at 1:3500 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit .

Exposure time: 180s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-03679-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-33421PEP	GABA-A R delta Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-03679

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

