

Product Datasheet

HDAC6 Antibody (236) - BSA Free NBP3-05295

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-05295

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-05295



NBP3-05295

HDAC6 Antibody (236) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	236
Preservative	15mM Sodium Azide
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4

Product Description	
Description	Novus Biologicals Mouse HDAC6 Antibody (236) - BSA Free (NBP3-05295) is a monoclonal antibody validated for use in WB, ELISA, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10013
Gene Symbol	HDAC6
Species	Human, Mouse (Negative), Zebrafish (Negative)
Specificity/Sensitivity	The mouse monoclonal antibody 236 recognizes D1 domain of human histone deacetylase 6 (HDAC6; intracellular antigen), amino acids 71-91. Crossreactivity with other species was not tested.
Immunogen	Human HDAC6


Product Application Details	
Applications	Western Blot, ELISA, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Flow (Intracellular)

Images


Western Blot: HDAC6 Antibody (236) [NBP3-05295] - Western blotting analysis of human HDAC6 by mouse monoclonal antibody NBP3-05295 in HEK-293T cell line under reducing conditions. PVDF membrane was blocked in 5% milk, and incubated with primary antibody (1 ug/ml), followed by HRP-conjugated anti-mouse secondary antibody.




Immunocytochemistry/Immunofluorescence: HDAC6 Antibody (236) [NBP3-05295] - Immunocytochemistry staining of HDAC6 in formaldehyde-fixed and Triton-permeabilized HEK-293T cells (A) and SH-SY5Y cells (B) by mouse monoclonal antibody NBP3-05295, followed by anti-mouse Alexa Fluor 488 (green), DNA indicated by DAPI (blue). In SH-SY5Y cells tubulin was stained by a rabbit polyclonal antibody, followed by anti-rabbit Alexa Fluor 647 (red, colocalization yellow).




ELISA: HDAC6 Antibody (236) [NBP3-05295] - ELISA analysis of human HDAC6 by mouse monoclonal antibody NBP3-05295. Plate was coated with recombinant human HDAC6, then blocked with BSA, and exposed to dilution series of anti-HDAC6 primary antibody, followed by HRP-conjugated anti-mouse secondary antibody, and colorimetric signal of processed OPD substrate was measured at 492 nm.



Immunoprecipitation: HDAC6 Antibody (236) [NBP3-05295] - Immunoprecipitation of recombinant human HDAC6 by mouse monoclonal antibodies 178 and NBP3-05295 using protein G-coated Dynabeads. Eluted proteins were processed by SDS PAAGE under reducing conditions and visualized by Coomassie Brilliant Blue.



Flow Cytometry: Mouse Monoclonal HDAC6 Antibody (236) [NBP3-05295] - Separation of K562 cells stained using anti-human HDAC6 (236) purified antibody (concentration in sample 5.33 ug/ml, GAM APC, red-filled) from K562 cells unstained by primary antibody (GAM APC, black-dashed) in flow cytometry analysis (intracellular staining).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-05295

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NB100-56343PEP	HDAC6 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-05295

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

