

# Product Datasheet

## Fibrillarin Antibody [DyLight 755] NBP3-05533IR

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-05533IR](http://www.novusbio.com/NBP3-05533IR)

Updated 12/13/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-05533IR](http://www.novusbio.com/reviews/destination/NBP3-05533IR)



**NBP3-05533IR**

Fibrillarin Antibody [DyLight 755]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Monoclonal  |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG1  |
| Conjugate                   | DyLight 755   |
| Purity                      | Immunogen affinity purified   |
| Buffer                      | 50mM Sodium Borate  |
| Product Description         |   |
| Host                        | Mouse   |
| Gene ID                     | 2091  |
| Gene Symbol                 | FBL   |
| Species                     | Human, Mouse, Rat   |
| Immunogen                   | Recombinant full length human Fibrillarin sequence expressed in and purified from E. coli.  |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.           |
| Product Application Details |   |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry                  |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-05533IR**

---

|               |  |
|---------------|--|
| NBP1-97005IR  | Mouse IgG1 Isotype Control (MG1) [DyLight 755] |
| NBP2-58592PEP | Fibrillarin Recombinant Protein Antigen        |
| D6050         | IL-6 [HRP]                                     |
| NBL1-10607    | Fibrillarin Overexpression Lysate              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-05533IR](http://www.novusbio.com/reviews/submit/NBP3-05533IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

