

# Product Datasheet

## LMAN2 Antibody (04) [FITC]

### NBP3-05831F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-05831F](http://www.novusbio.com/NBP3-05831F)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-05831F](http://www.novusbio.com/reviews/destination/NBP3-05831F)



**NBP3-05831F**

LMAN2 Antibody (04) [FITC]

| Product Information         |                                                                                                                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                       |
| Clonality                   | Monoclonal                                                                                                                                                                                                     |
| Clone                       | 04                                                                                                                                                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                             |
| Isotype                     | IgG                                                                                                                                                                                                            |
| Conjugate                   | FITC                                                                                                                                                                                                           |
| Purity                      | Protein A purified                                                                                                                                                                                             |
| Buffer                      | PBS                                                                                                                                                                                                            |
| Product Description         |                                                                                                                                                                                                                |
| Host                        | Mouse                                                                                                                                                                                                          |
| Gene ID                     | 10960                                                                                                                                                                                                          |
| Gene Symbol                 | LMAN2                                                                                                                                                                                                          |
| Species                     | Human                                                                                                                                                                                                          |
| Immunogen                   | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human LMAN2 (Uniprot#: Q12907; Asp45-Arg322). |
| Product Application Details |                                                                                                                                                                                                                |
| Applications                | ELISA, Sandwich ELISA Capture                                                                                                                                                                                  |
| Recommended Dilutions       | ELISA, Sandwich ELISA Capture                                                                                                                                                                                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-05831F**

---

|                    |                                              |
|--------------------|----------------------------------------------|
| NBP1-97278         | Mouse IgG Fc Isotype Control [FITC]          |
| H00010960-P01-10ug | Recombinant Human LMAN2 GST (N-Term) Protein |
| NBP2-07575         | LMAN2 Overexpression Lysate                  |
| NBP1-40256-0.1mg   | ATF6 Antibody (70B1413.1) - BSA Free         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-05831F](http://www.novusbio.com/reviews/submit/NBP3-05831F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

