

# Product Datasheet

## KLHL24 Antibody - BSA Free

### NBP3-05979-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-05979](http://www.novusbio.com/NBP3-05979)

Updated 2/26/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-05979](http://www.novusbio.com/reviews/destination/NBP3-05979)



**NBP3-05979-100ul**

KLHL24 Antibody - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.03% Proclin 300                                                                           |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Antigen and protein A Affinity-purified                                                     |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**


|                    |                                                                                                                                                                                                                             |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles. |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                      |
| <b>Gene ID</b>     | 54800                                                                                                                                                                                                                       |
| <b>Gene Symbol</b> | KLHL24                                                                                                                                                                                                                      |
| <b>Species</b>     | Human                                                                                                                                                                                                                       |
| <b>Immunogen</b>   | Produced in rabbits immunized with E. coli-derived Human KLHL24 fragment.                                                                                                                                                   |

**Product Application Details**

|                              |                                          |
|------------------------------|------------------------------------------|
| <b>Applications</b>          | Immunohistochemistry-Paraffin            |
| <b>Recommended Dilutions</b> | Immunohistochemistry-Paraffin 1:50-1:200 |

**Images**

Immunohistochemistry-Paraffin: KLHL24 Antibody [NBP3-05979] - Staining of human KLHL24 in human kidney with rabbit polyclonal antibody at 1:100 dilution.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-05979-100ul**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-37705PEP | KLHL24 Recombinant Protein Antigen                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-05979](http://www.novusbio.com/reviews/submit/NBP3-05979)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

