

# Product Datasheet

## Serine racemase Antibody [DyLight 405] NBP3-06588V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06588V](http://www.novusbio.com/NBP3-06588V)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06588V](http://www.novusbio.com/reviews/destination/NBP3-06588V)



**NBP3-06588V**

Serine racemase Antibody [DyLight 405]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Conjugate</b>     | DyLight 405                                                                                 |
| <b>Purity</b>        | Antigen and protein A Affinity-purified                                                     |
| <b>Buffer</b>        | 50mM Sodium Borate                                                                          |

**Product Description**

|                    |                                                                                    |
|--------------------|------------------------------------------------------------------------------------|
| <b>Host</b>        | Rabbit                                                                             |
| <b>Gene ID</b>     | 63826                                                                              |
| <b>Gene Symbol</b> | SRR                                                                                |
| <b>Species</b>     | Human                                                                              |
| <b>Immunogen</b>   | Produced in rabbits immunized with E. coli-derived Human Serine racemase fragment. |
| <b>Notes</b>       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Immunohistochemistry-Paraffin                                          |
| <b>Recommended Dilutions</b> | Immunohistochemistry-Paraffin                                          |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-06588V**

---

|                 |                                               |
|-----------------|-----------------------------------------------|
| NBP2-24891V     | Rabbit IgG Isotype Control [DyLight 405]      |
| NBP1-78868-10ug | Recombinant Human Serine racemase His Protein |
| NBL1-16460      | Serine racemase Overexpression Lysate         |
| NBP1-87069      | SDHB Antibody                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06588V](http://www.novusbio.com/reviews/submit/NBP3-06588V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

