

Product Datasheet

Recombinant Plant *Daucus carota* 1.0101 Strep (N-Term) Protein NBP3-06720

Unit Size: 0.1 mg

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-06720

Updated 5/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-06720



NBP3-06720Recombinant Plant *Daucus carota* 1.0101 Strep (N-Term) Protein


Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20 to -70C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Reconstitution Instructions	Reconstitute with sterile deionized water with or without azide.
Purity	>95%, by HPLC
Buffer	100 mM Tris, 150 mM NaCl, 0.01% Tween 20, pH 8.0
Target Molecular Weight	17.5 kDa

Product Description	
Description	<p>A recombinant protein with a N-Terminal Strep-tag corresponding to 167 amino acids of <i>Daucus carota</i> 1.0101</p> <p>Source: <i>E.coli</i></p> <p>Store product aliquots at -20C to -70C. Use stored aliquots within 6 months.</p>
Species	Plant
Reactivity Notes	Carrot (<i>Daucus carota</i>)
Preparation Method	The protein was purified by ionex and affinity chromatography, using Strep-tag. Endotoxin was removed using a specific endotrap carrier. Product was lyophilized after purification.


Product Application Details	
Applications	ELISA, Flow Cytometry
Recommended Dilutions	Flow Cytometry, ELISA

Images


Flow Cytometry: Recombinant Plant *Daucus carota* 1.0101 Strep (N-Term) Protein [NBP3-06720] - Flow cytometry dot-plot staining pattern of *Daucus carota* 1.0101 recombinant allergen stimulated human peripheral whole blood lymphocytes and basophils of a proven allergic donor stained using anti-human CD63 (NB100-77913) FITC and anti-human CD203c PE antibodies.



ELISA: Recombinant Plant *Daucus carota* 1.0101 Strep (N-Term) Protein [NBP3-06720] - ELISA (enzyme-linked immunosorbent assay) test was designed to prove the bond between the coated target recombinant allergen *Daucus carota* 1.0101 and allergen-specific human plasma IgG4 antibodies of *Daucus carota* positive donor. A measurable signal was subsequently generated by the addition of biotin labeled anti-human IgG4 (detection) antibody, Streptavidin-HRP and substrate solution (TMB). The intensity of the signal is proportional to the amount of coated *Daucus carota* 1.0101.



SDS-Page: Recombinant Plant *Daucus carota* 1.0101 Strep (N-Term) Protein [NBP3-06720] - Recombinant allergen *Daucus carota* 1.0101 purity verification. 5 µg of *Daucus carota* 1.0101 with > 95 % purity checked by Coomassie Brilliant Blue stained SDS-PAGE.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-06720

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

