

Product Datasheet

Mouse Calretinin ELISA Kit (Colorimetric) NBP3-06813

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-06813

Updated 9/23/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-06813



NBP3-06813**Mouse Calretinin ELISA Kit (Colorimetric)**


| Product Information | |
|----------------------------|--|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Storage of components varies. See protocol for specific instructions. |

| Product Description | |
|--------------------------------|---|
| Gene ID | 794 |
| Gene Symbol | CALB2 |
| Species | Mouse |
| Specificity/Sensitivity | This kit recognizes Mouse Calretinin in samples. No significant Calretininoss-reactivity or interference between Mouse Calretinin and analogues was observed. |
| Kit Components | MiCalretinino ELISA Plate(Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product DesCalretininipion, Certificate of Analysis |
| Standard Curve Range | 0.16 - 10 ng/mL |
| Sensitivity | 0.09 ng/mL |
| Inter Assay | 4.74 |
| Intra Assay | 3.85 |
| Assay Type | Sandwich-ELISA |
| Spiking Recovery | 88 - 111 % |
| Suitable Sample Type | Serum, plasma, cell supernatant, and other biological fluids. |
| Sample Volume | 100 ul |

| Product Application Details | |
|------------------------------------|-----------------------|
| Applications | ELISA, Sandwich ELISA |
| Recommended Dilutions | ELISA, Sandwich ELISA |

Images

ELISA: Mouse Calretinin ELISA Kit (Colorimetric) [NBP3-06813] - Standard Curve Reference.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-06813

| | |
|------------------|--|
| NBP1-48344-0.1mg | Recombinant Human Calretinin His Protein |
| MAB1455 | Albumin Antibody (188835) [Unconjugated] - Serum |
| NBP1-88220 | Calretinin Antibody |
| NB120-11427 | Parvalbumin Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-06813

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

