

# Product Datasheet

## Rat CX3CR1 ELISA Kit (Colorimetric) NBP3-06833

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06833](http://www.novusbio.com/NBP3-06833)

Updated 9/23/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06833](http://www.novusbio.com/reviews/destination/NBP3-06833)



**NBP3-06833****Rat CX3CR1 ELISA Kit (Colorimetric)**


| <b>Product Information</b> |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 1 Kit                                                                                                    |
| <b>Concentration</b>       | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| <b>Storage</b>             | Storage of components varies. See protocol for specific instructions.                                    |

| <b>Product Description</b>     |                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>                 | 1524                                                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Symbol</b>             | CX3CR1                                                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>                 | Rat                                                                                                                                                                                                                                                                                                                                                                |
| <b>Specificity/Sensitivity</b> | This kit recognizes Rat CX3CR1 in samples. No significant cross-reactivity or interference between Rat CX3CR1 and analogues was observed.                                                                                                                                                                                                                          |
| <b>Kit Components</b>          | Micro ELISA Plate(Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Description, Certificate of Analysis |
| <b>Standard Curve Range</b>    | 1.56 - 100 ng/mL                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sensitivity</b>             | 0.94 ng/mL                                                                                                                                                                                                                                                                                                                                                         |
| <b>Inter Assay</b>             | 5.34                                                                                                                                                                                                                                                                                                                                                               |
| <b>Intra Assay</b>             | 5.45                                                                                                                                                                                                                                                                                                                                                               |
| <b>Assay Type</b>              | Sandwich-ELISA                                                                                                                                                                                                                                                                                                                                                     |
| <b>Spiking Recovery</b>        | 83 - 108 %                                                                                                                                                                                                                                                                                                                                                         |
| <b>Suitable Sample Type</b>    | Serum, plasma, cell supernatant, and other biological fluids.                                                                                                                                                                                                                                                                                                      |
| <b>Sample Volume</b>           | 100 ul                                                                                                                                                                                                                                                                                                                                                             |

| <b>Product Application Details</b> |                       |
|------------------------------------|-----------------------|
| <b>Applications</b>                | ELISA, Sandwich ELISA |
| <b>Recommended Dilutions</b>       | ELISA, Sandwich ELISA |

**Images**

ELISA: Rat CX3CR1 ELISA Kit (Colorimetric) [NBP3-06833] - Standard Curve Reference.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-06833**

---

|                     |                                  |
|---------------------|----------------------------------|
| NBP1-76949PEP-0.1mg | CX3CR1 Antibody Blocking Peptide |
| 210-TA-005          | TNF-alpha [Unconjugated]         |
| AF5825              | CX3CR1 Antibody [Unconjugated]   |
| M6000B-1            | IL-6 [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06833](http://www.novusbio.com/reviews/submit/NBP3-06833)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

