

# Product Datasheet

## **p16INK4a/CDKN2A Antibody [DyLight 755] NBP3-06986IR**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06986IR](http://www.novusbio.com/NBP3-06986IR)

Updated 7/11/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06986IR](http://www.novusbio.com/reviews/destination/NBP3-06986IR)



**NBP3-06986IR**

p16INK4a/CDKN2A Antibody [DyLight 755]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 755
Purity	Peptide affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	1029
Gene Symbol	CDKN2A
Species	Human
Reactivity Notes	[Species Reactivity]
Specificity/Sensitivity	Human CDKN2A has 6 isoforms, including isoform 1 (156aa, 17kD), isoform 2 (105aa, 11kD), isoform 3 (116aa, 12kD), isoform 5 (167aa, 18kD), isoform tumor suppressor ARF (132aa, 14kD) and isoform smARF (85aa, 9kD). Mouse CDKN2A has 4 isoforms, including isoform 1 (168aa, 18kD), isoform 2 (125aa, 14kD), isoform tumor suppressor ARF (169aa, 19kD) and isoform smARF (125aa, 14kD). Rat CDKN2A has two isoforms, including isoform 1 (159aa, 17kD) and isoform ARF (160aa, 17kD). NBP3-06986 can detect only human isoforms 1 and 5, but cannot detect mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 17 amino acids near the carboxy terminus of human CDKN2A. The immunogen is located within the last 50 amino acids of CDKN2A.
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-06986IR**

---

NBP2-24891IR	Rabbit IgG Isotype Control [DyLight 755]
NBP2-35199-10ug	Recombinant Human p16INK4a/CDKN2A Protein
NBP2-35200-5ug	Recombinant Human p16INK4a/CDKN2A TAT Protein
NB200-103	p53 Antibody (PAb 240) - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06986IR](http://www.novusbio.com/reviews/submit/NBP3-06986IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

