

Product Datasheet

p16INK4a/CDKN2A Antibody [mFluor Violet 500 SE] NBP3-06986MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-06986MFV500

Updated 10/24/2023 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-06986MFV500



NBP3-06986MFV500

p16INK4a/CDKN2A Antibody [mFluor Violet 500 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	mFluor Violet 500 SE
Purity	Peptide affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	1029
Gene Symbol	CDKN2A
Species	Human
Specificity/Sensitivity	Human CDKN2A has 6 isoforms, including isoform 1 (156aa, 17kD), isoform 2 (105aa, 11kD), isoform 3 (116aa, 12kD), isoform 5 (167aa, 18kD), isoform tumor suppressor ARF (132aa, 14kD) and isoform smARF (85aa, 9kD). Mouse CDKN2A has 4 isoforms, including isoform 1 (168aa, 18kD), isoform 2 (125aa, 14kD), isoform tumor suppressor ARF (169aa, 19kD) and isoform smARF (125aa, 14kD). Rat CDKN2A has two isoforms, including isoform 1 (159aa, 17kD) and isoform ARF (160aa, 17kD). NBP3-06986 can detect only human isoforms 1 and 5, but cannot detect mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 17 amino acids near the carboxy terminus of human CDKN2A. The immunogen is located within the last 50 amino acids of CDKN2A.
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-06986MFV500

NBP2-35199-10ug	Recombinant Human p16INK4a/CDKN2A Protein
NBP2-35200-5ug	Recombinant Human p16INK4a/CDKN2A TAT Protein
NB200-103	p53 Antibody (PAb 240) - BSA Free
NB100-2736	MDM2/HDM2 Antibody (SMP14) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-06986MFV500

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

