

Product Datasheet

SARS-CoV-2 Nucleocapsid Antibody (3868) - BSA Free NBP3-07070

Unit Size: 0.05 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-07070

Updated 10/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-07070




NBP3-07070**SARS-CoV-2 Nucleocapsid Antibody (3868) - BSA Free**


Product Information	
Unit Size	0.05 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3868
Preservative	0.1% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	0.01M PBS (pH 7.2)
Product Description	
Description	This antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
Host	Mouse
Gene ID	43740575
Gene Symbol	N
Species	SARS-CoV-2
Specificity/Sensitivity	Specific for the nucleocapsid of Covid-19 virus. Cross reacts with SARS. These antibodies are non-reactive with: Influenza A, Influenza B, Adenovirus and MERS.
Immunogen	This antibody was developed using a recombinant, full-length nucleoprotein.
Product Application Details	
Applications	ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	ELISA 1:20-1:200, Immunocytochemistry/ Immunofluorescence 1:10-1:50
Application Notes	Sandwich ELISAs for SARS-CoV-2 Nucleocapsid Matched Pair Antibodies. A sandwich ELISA was performed using SARS-CoV-2 Nucleocapsid antibody (NBP3-07038, 2ug/ml) as capture antibody, a nucleocapsid recombinant protein as the binding protein, and the anti-SARS-CoV-2 Nucleocapsid antibody (NBP3-07070, 0.5ug/ml) as the detection antibody. Secondary: Goat anti-mouse IgG HRP conjugate at 1:20000 dilution. Detection range is from 0.03 ng to 300 ng. EC50 = 4.27 ng. A sandwich ELISA was performed using SARS-CoV-2 Nucleocapsid antibody (NBP3-07039, 2ug/ml) as capture antibody, the Nucleocapsid recombinant protein as the binding protein, and the anti-SARS-CoV-2 Nucleocapsid antibody (NBP3-07068, 0.5ug/ml) as the detection antibody. Secondary: Goat anti-mouse IgG HRP conjugate at 1:20000 dilution. Detection range is from 0.03 ng to 300 ng. EC50 = 15.82 ng.

Images

Sandwich ELISA: SARS-CoV-2 Nucleocapsid Antibody (3868) [NBP3-07070] - Sandwich ELISA for SARS-CoV-2 Nucleocapsid Matched Pair Antibodies. A sandwich ELISA was performed using SARS-CoV-2 Nucleocapsid antibody (NBP3-07038, 2ug/ml) as capture antibody, a nucleocapsid recombinant protein as the binding protein, and the anti-SARS-CoV-2 Nucleocapsid antibody (NBP3-07070, 0.5ug/ml) as the detection antibody. Secondary: Goat anti-mouse IgG HRP conjugate at 1:20000 dilution. Detection range is from 0.03 ng to 300 ng. EC50 = 4.27 ng



Sandwich ELISA: SARS-CoV-2 Nucleocapsid Antibody (3868) [NBP3-07070] - Sandwich ELISA for SARS-CoV-2 Nucleocapsid Antibodies. A sandwich ELISA was performed using SARS-CoV-2 Nucleocapsid antibody (NBP3-07039, 2ug/ml) as capture antibody, the Nucleocapsid recombinant protein as the binding protein, and the anti-SARS-CoV-2 Nucleocapsid antibody (NBP3-07068, 0.5ug/ml) as the detection antibody. Secondary: Goat anti-mouse IgG HRP conjugate at 1:20000 dilution. Detection range is from 0.03 ng to 300 ng. EC50 = 15.82 ng



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-07070

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-90975	Recombinant SARS-CoV-2 Nucleocapsid His (C-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-07070

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

