

Product Datasheet

HR6B/UBE2B Antibody (PCRP-UBE2B-1C7) NBP3-07169-100ug

Unit Size: 100 ug

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-07169

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-07169



NBP3-07169-100ug

HR6B/UBE2B Antibody (PCRP-UBE2B-1C7)

Product Information	
Unit Size	100 ug
Concentration	0.2 mg/ml
Storage	Store at 4C.
Clonality	Monoclonal
Clone	PCRP-UBE2B-1C7
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Purity	Protein A or G purified
Buffer	10 mM PBS with 0.05% BSA
Target Molecular Weight	17 kDa


Product Description	
Description	<p>200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-08912)</p> <p>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C.</p>
Host	Mouse
Gene ID	7320
Gene Symbol	UBE2B
Species	Human
Reactivity Notes	Predicted to work in Mouse, Rat, Xenopus and Zebrafish.
Immunogen	Recombinant full-length human HR6B/UBE2B protein (Uniprot: P63146)

Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Protein Array
Recommended Dilutions	Flow Cytometry 1-2 ug/million cells, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Protein Array
Application Notes	Optimal dilution for a specific application should be determined.




Images

Immunocytochemistry/Immunofluorescence: HR6B/UBE2B Antibody (PCRP-UBE2B-1C7) [NBP3-07169] - Immunofluorescent Analysis of PFA-fixed HeLa cells. HR6B/UBE2B Mouse Monoclonal Antibody (PCRP-HR6B/UBE2B-1C7) followed by IgG-CF488 (green), counterstained with phalloidin.



Flow Cytometry: HR6B/UBE2B Antibody (PCRP-UBE2B-1C7) [NBP3-07169] - Flow Cytometric Analysis of PFA-fixed HeLa cells. HR6B/UBE2B Mouse Monoclonal Antibody (PCRP-HR6B/UBE2B-1C7) followed by goat anti- Mouse IgG-CF488 (blue); unstained cells (red).



Protein Array: HR6B/UBE2B Antibody (PCRP-UBE2B-1C7) [NBP3-07169] - Analysis of Protein Array containing more than 19,000 full-length human proteins using HR6B/UBE2B Mouse Monoclonal Antibody (PCRP-HR6B/UBE2B-1C7).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-07169-100ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NBP2-35088-10ug	Recombinant Human HR6B/UBE2B His (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-07169

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

